



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Isolation Station, Wall or Door Mount, Standard Package #ISO2825

Specifications



Construction:

- Painted aluminum

Features:

- Saves floor space, no more carts getting in the way of patient flow
- Allows storage of PPE (personal protective equipment) in areas with space limitations
- Pre-drilled holes for wall mount (mounting hardware included)
- Door-hanging bracket included to hang on 2" door
- Designed to hold 3 boxes of gloves, bulk pack of disposable gowns, face masks or shoe covers and caps
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish

Other Wall Mount / Door Hanging Isolation Station Options:

- Silver ion antimicrobial paint upgrade available
- Custom configurations are available based on your requirements

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Overall Dimensions:

- Wall mount: 22" H x 24.9" W x 3.7"D (55.9 cm H x 63.2 cm W x 9.4 cm D)
- Net weight: 15 pounds (6.8 kg)